

## 2 ml Screw Vials with 9 mm Screw Cap

### Agilent Technologies compatible (HP) 9 mm Screw Cap

The 9 mm screw thread was especially developed for the Agilent Autosampler to allow the handling with the robotic arm. When using screw-top Vials, we generally recommend the Agilent compatible screw cap with 9 mm thread as they seal best from all the available screw-tops. To our knowledge, they are compatible with all modern autosamplers.



#### G004-HP-H G004-HP-D

2 ml Agilent compatibles Screw Vial,  
clear (H) or amber (D)  
borosilicate glass;  
size : 12 x 32 mm

#### 9 mm Screw Caps

##### with ms-Pure non-pigmented Silicone/PTFE Septum

The ms-Pure Silicone/PTFE Septum is made from non-pigmented Silicone/PTFE. Tests show that the MS-Pure Septum is over 60 % cleaner than the standard white/red versions on the market. With the increasing sensitivity of GC's and HPLC's there is an increased need for low bleed septum material to reduce the potential for sample error and to improve overall reliability.



### Our autumn offer valid till 31.10.2011



#### G001X-HP-SKFK10-H

clear glass: combination pack  
2 x 500 pcs. G004-HP-H  
2 x 500 pcs. G004-HP-C\*-SKFK10



**combination price CHF 28.50/100 pcs.**

*instead of CHF 34.50/100 pcs.*

#### G001X-HP-SKFK10-D

amber glass: combination pack  
2 x 500 pcs. G004-HP-D  
2 x 500 pcs. G004-HP-C\*-SKFK10



**combination price CHF 28.90/100 pcs.**

*instead of CHF 35.00/100 pcs.*

Cap colour depending on stock  
offer valid as long as stock lasts; a maximum of 25.000 pcs./customer  
Price excl. VAT, excl. Postage and Packaging CHF 17.50